Ionoscell

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Ionoscell

Glass ionomer permanent luting cement

Properties

IONOSCELL is a fast high refractive index glass. It contains vacuum dried polyacrylic acid. When mixed with water it produces a luting cement for the permanent scaling of dental prostheses.

Directions for use

Since the powder tends to settle, it is recommended to shake the bottle so as to fluff powder. Fill the measuring spoon to level with powder. The water dispensing bottle should be half filled with distilled water and held vertically downwards to allow dripping.

The powder should be placed separately on the mixing pad. A first portion of powder should be incorporated into water. As soon as it is completely wet, add the second portion and mix quickly so as to obtain a ready for use smooth cream. The cement in this glossy state should be applied immediately to the prosthetic element to be sealed. Then seat the prosthesis on the dry prepared tooth. The excess cement should be carefully removed when still rubbery before complete curing.

Mixing ratio

2 measuring spoonful + 3 drops of water

Mixing time : 15 seconds

Working time : 2 minutes

Setting time : 2 minutes

Total time : 4 to 5 minutes at 23°C

Precautions for use

This is an adhesive cement. Mixing pad and spatula should be wiped clean with damp napkin just prior to setting.

The cement contains PAA. Avoid contact with eyes. In case of eye contact, rinse immediately with plenty of water and seek medical advice.

Harmful if swallowed.

Store in a dry place. The lid must be tightly closed immediately after use.

Presentation

Kit containing :

- 30 g powder,
- water dispensing bottle,
- mixing pad,
- measuring spoon.

For professional use only in the practice of dentistry.

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