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EN

PRODUCT DESCRIPTION

RESORCELL Gauze Absorbable Oxidized Cellulose Gauze Hemostat is a sterile, biodegradable, knitted gauze, prepared by the controlled oxidation of cellulose. It is intended for use as topical hemostat and serves as a hemostatic adjunct in the control of hemorrhage. The gauze is white and can be sutured or cut without fraying. In use RESORCELL Gauze swells minimally as blood and other fluids are absorbed.

The bacteriostatic / bactericidal properties of RESORCELL Gauze inhibit aerobic and anaerobic bacteria cell growth and multiplication.

Although the product has demonstrated bacteriostatic / bactericidal properties, it is not intended nor claimed as a substitute for systemically administered therapeutic or prophylactic anti-microbial agent to control or prevent post-operative infections.

It is packaged individually in (double) blister, sterilized by means of Gamma Irradiation, and is for single use only.

INDICATIONS FOR USE

Topical hemostat for use as an adjunct to hemostasis by tamponade effect after tooth extraction, in particular where control of capillary, venous, and arteriolar bleeding, by pressure, ligature, and other conventional procedures, is either ineffective or impractical.

CONTRAINDICATIONS

RESORCELL Gauze should not be used for implantation in bone defects, such as fractures, since there is a possibility of interference with callus formation and a theoretical chance of cyst formation. RESORCELL Gauze should not be used in bleeding from large arteries.

RESORCELL Gauze should not be used in chemically cauterized areas; its use should not be preceded by application of escharotic chemicals (silver nitrate or any other).

RESORCELL Gauze should not be used as a surface dressing except for immediate control of bleeding as it inhibits new skin growth.

RESORCELL Gauze should not be used in closure of skin incisions as it may interfere with healing of the skin edges. This is due to mechanical interposition of oxidized cellulose and not to intrinsic interference with wound healing.

SIDE EFFECTS

Occasional reports of "burning" and "stinging" sensations are believed to be due to the low pH of the product.

DIRECTIONS FOR USE

Sterile technique should be applied.

RESORCELL Gauze should be applied dry and should not be moistened prior to application.

RESORCELL Gauze is best employed by loose application against the bleeding surface, where it should adhere.

Remove any excess prior to surgical closure to facilitate bio-degradation and to minimize the possibility of foreign body reaction! Bio-degradation depends upon several factors including the amount used, degree of saturation with blood and the tissue bed. RESORCELL Gauze can be cut to size.

IMPORTANT NOTICE !

Surgeons should be aware that encapsulation of fluid and foreign body reactions have been reported for this type of product. During follow up screening or imaging this should be considered as a possible cause of certain clinical findings that may be otherwise suggestive of pathology.

RESORCELL Gauze expands minimally with the absorption of

fluids. In closed tissue spaces with presence of bone, RESORCELL Gauze should be removed when bleeding has stopped, if possible. Failure to do so might result in unintended pressure on neighbouring structures which may cause pain for the patient or could create the potential for nerve damage.

While most reports have been in connection with laminectomy, reports of paralysis have also been received in connection with other procedures.

By absorbing fluid RESORCELL Gauze may expand and impinge on neighbouring structures. Therefore care should be taken to avoid over-packing, particularly within bony cavities, as it may lead to interference with normal function and may eventually lead to compression necrosis of surrounding tissue.

In cases of postoperative infection, re-operation may be necessary to remove infected material and allow drainage.

WARNINGS

RESORCELL Gauze is not intended as a substitute for careful surgery and the proper use of sutures and ligatures. Closing the product in a contaminated wound without drainage may lead to complications and should be avoided.

Avoid use in infected areas.

There have been reports of stenotic effect when oxidized cellulose has been applied as a wrap during vascular surgery. Although it has not been established that the stenosis was directly related to the use of oxidized cellulose, it is important to be cautious and avoid applying the material lightly as a wrapping.

Care should be taken to avoid the introduction of fragments of the material directly into the circulatory system.

RESORCELL Gauze should not be impregnated with anti-infective agents or with other materials such as buffering or hemostatic substances. Its hemostatic effect is not enhanced by the addition of thrombin, the activity of which is destroyed by the low pH of the product.

Special care must be taken by practitioner, regardless of the type of surgical procedure, to consider the advisability of removing RESORCELL Gauze after hemostasis is achieved.

GENERAL STORAGE AND HANDLING

RESORCELL Gauze should be stored dry in its sealed protective outer packaging. Keep away from direct sunlight. For storage temperature conditions please refer to label on box and pouches. Warning! Browning of the gauze and tendency to crumble indicate incorrect storage.

Remove the product from the sterile package using sterile technique. Prior to opening, the integrity of the sterile barrier should be confirmed.

Do not use if packaging has been damaged.

Any unused contents should be discarded.

The packaging is subject to and conforms to European & International legislation and other applicable standards. The packaging protects the product from external influences and guarantees its sterility during storage and transport, under normal use.

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